



MIGRATING ELN ENTRIES FROM TEST TO THE PRODUCTION SERVER

Dr. Manuel Ehling Scientific Infrastructure MDC 15.12.2022






Go to Workspace on the TEST server and choose the folders you want to migrate




Click on Export








Choose the .ZIP bundle containing .XML files





Recommendation: Choose depth of infinity to export Click on Export to finalize




HELMHOLTZ

12/14/2022




Wait until you receive a message from the system (see bell sign; takes a few seconds) and download ZIP-file from the provided link



COMMENT




There is also the option to Export all of your files from test to production However, this would also include the test documents you created in the past





Go to the PRODUCTION server or any other receiving version of RSpace




Please note:

There are no folders of mine on this server, yet




Go to My RSpace and the Export/Import window on the right side





Choose the downloaded files to re-upload and then press Import





Check for error messages or other potential issues with your data files



Workspace Gallery Messaging Apps Inventory My RSpace Account ▼

Archive Import Report

Back to workspace Import another archive

Archive import was successful.

Validation results

Rule	Result	Description
ZIP_FILE	PASS	The archive should be a zip file
CHECKSUM_CALCULATED	PASS	RSpace needs to be able generate a checksum of the archive.
MANIFEST_FILE_PRESENT	PASS	The archive must contain a manifest file called manifest.txt
MANIFEST_SOURCE	UNTESTED	The archive manifest must contain a property 'archiveSource'
FOLDER_TREE_PRESENT	PASS	The archive must contain an XML file of the folder tree called 'folderTree.xml'
CHECKSUM_MATCHES	UNTESTED	The checksum of the archive must equal the checksum at the time of archive creation.
XMLSCHEMA	PASS	XML files must conform to their XML schemas
GENERAL_ARCHIVE_STRUCTURE	UNTESTED	The archive should contain all required resources.
USER_SCHEMA_VERSION_RANGE_OK	UNTESTED	The User schema is incompatible with the current database version
USER_FILE_READABLE	UNTESTED	The users.xml file could not be parsed
DOC_SCHEMA_VERSION_RANGE_OK	PASS	The Document schema of at least one record is incompatible with the current database version
FORM_SCHEMA_VERSION_RANGE_OK	UNTESTED	The Form schema of at least one form is incompatible with the current database version
FOLDER_SCHEMA_VERSION_RANGE_OF	(PASS	The Folder schema of at least one form is incompatible with the current database version
FOLDER_FILE_READABLE	PASS	The folders.xml file could not be parsed
UNKNOWN	UNTESTED	The archive cannot have unexpected i did content
ARCHIVE_NOT_TOO_NEW	UNTESTED	The archive was made from a newer version of RSpace than this one

Documents imported


Name	Global I
Project description and aim	SD1959
Untitled document	SD19592
Cell proliferation as ray	SD19593
Errors	

HELMHOLTZ 11 12/14/2022 Info



12

Now the imported data files are visible under my workspace




Please note:

Now my newly imported folder "Project C" is visible under my workspace



Please recheck briefly whether all files were re-imported correctly




ADDITIONAL COMMENTS



Linked images will be added to your new gallery directly

However, it is also possible to import images and folders from the gallery independently



ADDITIONAL COMMENTS



- Linked images will be added to your new gallery directly
- However, it is also possible to import images and folders from the gallery independently
- Documents, which were shared on the test installation have to be shared again
- Settings of the lab have to be setup again (if not done, yet)

ARE THERE OPEN QUESTIONS OR ADDITIONAL FEEDBACK?

